

Publications (Apr.2012 - Mar. 2013)

1. R. Asaoka, T. Uemura, S. Nishida, T. Fujiwara, T. Ueda and A. Nakano: "New insights into the role of Arabisopsis RABA1 GTPases in salinity stress tolerance," *Plant Signal. Behav.* in press.
2. T. Jia, C. Gao, Y. Cui, J. Wang, Y. Ding, Y. Cai, T. Ueda, A. Nakano and L. Jiang: "ARA6 (Q69L) expression in transgenic Arabidopsis cells induces the formation of enlarged multivesicular bodies," *J. Exp. Biol.* in press.
3. Y. Ueda, A. Makino, K. Murase-Tamada, S. Sakai, T. Inaba, F. Hullin-Matsuda, and T. Kobayashi: "Sphingomyelin regulates the transbilayer movement of diacylglycerol in the plasma membrane of Madin-Darby canine kidney cells", *FASEB J.* in press.
4. S.-W. Choi, T. Tamaki, K. Ebine, T. Uemura, T. Ueda, and A. Nakano: "Distinct RABA members act in multiple steps of subcellular trafficking of the FLAGELLIN SENSING 2 receptor," *Plant Cell* 25, 1174-1187 (2013).
5. T. Uejima, K. Ihara, M. Sunada, M. Kawasaki, T. Ueda, R. Kato, A. Nakano, and S. Wakatsuki: "Direct metal recognition by guanine nucleotide exchange factor in the initial step of exchange reaction," *Acta Cryst. D* 69, 345-351 (2013).
6. R. Asaoka, T. Uemura, J. Ito, M. Fujimoto, E. Ito, T. Ueda, and A. Nakano: "Arabidopsis RABA1 GTPases are involved in transport between the *trans*-Golgi network and the plasma membrane, and are required for salinity stress tolerance," *Plant J.* 73, 240-249 (2013).
7. A. Era, N. Kutsuna, T. Higaki, S. Hasezawa, A. Nakano, and T. Ueda: "Microtubule stability affects the unique motility of F-actin in *Marchantia polymorpha*," *J. Plant Res.* 126, 113-119 (2013).
8. K. Ebine, T. Uemura, A. Nakano, and T. Ueda: "Flowering time modulation by a vacuolar SNARE via *FLOWERING LOCUS C* in *Arabidopsis thaliana*," *PLoS ONE*. 7, e442239 (2012).
9. T. Uemura, T. Ueda, and A. Nakano: "The physiological role of SYP4 in the salinity and osmotic stress tolerances," *Plant Signal. Behav.* 7, 1118-1120 (2012).
10. Y. Ito, T. Uemura, K. Shoda, M. Fujimoto, T. Ueda, and A. Nakano: "*cis*-Golgi proteins accumulate near the ER exit sites and act as the scaffold for Golgi regeneration after brefeldin A treatment in tobacco BY-2 cells," *Mol. Biol. Cell* 23, 3203-3214 (2012).
11. E. Feraru, M. I. Feraru, R. Asaoka, T. Paciorek, R. De Rycke, H. Tanaka, A. Nakano, and J. Friml: "BEX5/RabA1b acts as a regulator of the TGN-to-PM protein trafficking in *Arabidopsis*," *Plant Cell* 24, 3074-3086 (2012).
12. T. Hamada, M. Tominaga, T. Fukaya, M. Nakamura, A. Nakano, Y. Watanabe, T. Hashimoto, and T. I. Baskin: "RNA processing bodies, peroxisomes, Golgi bodies, mitochondria, and ER tubule junctions frequently pause at cortical microtubules," *Plant Cell Physiol.* 53, 699-708 (2012).
13. M. Okamoto, K. Kurokawa, K. Matsuura-Tokita, C. Saito, R. Hirata, and A. Nakano: "High-curvature domains of the endoplasmic reticulum (ER) are important for the organization of ER exit sites in *Saccharomyces cerevisiae*," *J. Cell Sci.* 125, 3412-3420 (2012).
14. K. Ebine, N. Miyakawa, M. Fujimoto, T. Uemura, A. Nakano, and T. Ueda: "Endosomal trafficking pathway regulated by ARA6, a RAB5 GTPase unique to plants," *Small GTPases* 3, 23-27 (2012).
15. S.-I. Morita, S. Takanezawa, A. Date, S. Watanabe, T. N. Murakami, N. Kawashima, H. Sato, and Y. Sako: "Raman spectroscopic analysis of H₂O₂-stimulated three-dimensional human skin models," *J. Spectroscopy* 2013, 903450 (6 pages) (2012).
16. H.-S. Shin, A. Okamoto, Y. Sako, S. Kim, C.-G. Pack, and S. Y. Kim: "Characterization of the triplet state of hybridization-sensitive DNA probe by using fluorescence correlation spectroscopy," *J. Phys. Chem. A.* 117, 27-33 (2013).
17. K. Oishi, H. Kurahashi, C.-G. Pack, Y. Sako, and Y. Nakamura: "A bipolar functionality of Q/N-rich proteins: Lsm4 amyloid causes clearance of yeast prions," *MicrobiologyOpen*, doi: 10.1002/mbo3.83 (2013).
18. N. Imamoto: "Cargo recognition explains nuclear transport regulation induced by nuclear pore complex reorganization," *J. Mol. Biol.* 425, 1849-2851 (2013).

19. T. Chinen, S. Kazami, Y. Nagumo, I. Hayakawa, A. Ikeda, M. Takagi, A. Yokosuka, N. Imamoto, Y. Mimaki, H. Kigoshi, H. Osada, and T. Usui: "Glaziovianin A prevents endosome maturation via inhibiting microtubule dynamics," *ACS Chem. Biol.* 8, 884-889 (2013).
20. M. Clever, Y. Mimura, T. Funakoshi, and N. Imamoto: "Regulation and coordination of nuclear envelope and nuclear pore complex assembly," *Nucleus* 4, 105-114 (2013).
21. M. Kimura, S. Kose, N. Okumura, K. Imai, M. Furuta, N. Sakiyama, K. Tomii, P. Horton, T. Takao, and N. Imamoto: "Identification of cargo proteins specific for the nucleocytoplasmic transport carrier transportin by combination of an in vitro transport system and stable isotope labeling by amino acids in cell culture (SILAC)-based quantitative proteomics," *Mol. Cell. Proteomics* 12, 145-157 (2013).
22. S. Hihara, C.G. Pack, K. Kaizu, T. Tani, T. Hanafusa, T. Nozak, S. Takemoto, T. Yoshimi, H. Yokota, N. Imamoto, N. Y. Sako, M. Kinjo, K. Takahashi, T. Nagai, and K. Maeshima: "Local nucleosome dynamics facilitate chromatin accessibility in living Mammalian cells," *Cell Rep.* 2, 1645-1656 (2012).
23. U. Kato, H. Inadome, M. Yamamoto, K. Emoto, T. Kobayashi, and M. Umeda: "Role for phospholipid flippase complex of ATP8A1 and CDC50A proteins in cell migration," *J. Biol. Chem.* 288: 4922 (2013).
24. M. Kudo, S. Wada, T. Ogawa, S. Namiki, A. Shinjo, K. Suzuki, Y. Ozawa, and K. Tsubota: "Confocal micro-imaging system for organic materials" The Third International Conference on nanoManufacturing nanoMan2012, 217-220, (2012).
25. R. Yachi, Y. Uchida, BH. Balakrishna, G. Anderluh, T. Kobayashi, T. Taguchi, and H. Arai: "Subcellular localization of sphingomyelin revealed by two toxin-based probes in mammalian cells", *Genes Cells* 17, 720-727 (2012).
26. J. M. Alam, T. Kobayashi, and M. Yamazaki: "Single giant unilamellar vesicle method reveals lysenin-induced pore formation in lipid membranes containing sphingomyelin," *Biochemistry* 51, 5160-5172 (2012).
27. F. Hullin-Matsuda, N. Tomishige, S. Sakai, R. Ishitsuka, K. Ishii, A. Makino, P. Greimel, M. Abe, E. L. Laviad, M. Lagarde, H. Vidal, T. Saito, H. Osada, K. Hanada, A. H. Futterman, and T. Kobayashi: "Limonoid compounds inhibit sphingomyelin biosynthesis by preventing CERT protein-dependent extraction of ceramides from the endoplasmic reticulum," *J. Biol. Chem.* 287, 24397-24411 (2012).
28. L. D. Chiu, F. Hullin-Matsuda, T. Kobayashi, H. Torii, H. O. Hamaguchi: "On the origin of the 1602 cm(-1) Raman band of yeasts; contribution of ergosterol," *J. Biophotonics* 1-5 DOI 10.1002/jbio.201200020 (2012).
29. J. Ikenouchi, M. Suzuki, K. Umeda, K. Ikeda, R. Taguchi, Te. Kobayashi, S. B. Sato, To. Kobayashi, D. B. Stoltz, and M. Umeda: "Lipid polarity is maintained in the absence of tight junctions," *J. Biol. Chem.* 287, 9525-9533 (2012).
30. M. Takahashi, T. Shibata, T. Yanagida, and Y. Sako: "A protein switch with tunable steepness reconstructed in *E. coli* cells with eukaryotic signaling proteins," *Biochem. Biophys. Res. Comm.*, 421, 731-735 (2012).
31. M. Ito, A. Shibata, J. Zhang, M. Hiroshima, Y. Sako, Y. Nakano, K. Kojima-Aikawa, B. Mannervik, S. Shuto, Y. Ito, R. Morgenstern, and H. Abe: Universal caging group for in cell detection of glutathione transferase applied to ¹⁹F NMR and bioluminogenic probes," *Chem. Bio. Chem.* 13, 1428-1432 (2012).
32. Y. Haga, K. Ishii, K. Hibino, Y. Sako, Y. Ito, N. Taniguchi, and T. Suzuki: "Visualizing specific protein glycoforms by transmembrane fluorescence resonance energy transfer," *Nat. Commun.* 3, 907 (2012).
33. H.-S. Shin, C.-G. Pack, S. Y. Kim, Y. Sako, and A. Okamoto: "Radiationless deactivation of doubly thiazole orange labeled DNA," *J. Luminescence* 130, 1784-1786 (2012).
34. M. Hiroshima, Y. Saeki, M. Okada-Hatakeyama, and Y. Sako: "Dynamically varying interactions between heregulin and ErbB proteins detected by single-molecule analysis in living cells," *Proc. Natl. Acad. Sci. USA* 109, 13984-13989 (2012).
35. C.-G. Pack, M.-R. Song, E.-J. Lee, M. Hiroshima, K.-H. Byun, J.-S. Kim, and Y. Sako: "Microenvironments and different nanoparticle dynamics in living cells revealed by a standard nanoparticle," *J. Cont. Rel.* 163(3), 315-321(2012).
36. K. Okamoto, and Y. Sako: "Variational Bayes analysis of a photon-based hidden Markov model for

single-molecule FRET trajectories,” *Biophys. J.* 103, 1315-1324 (2012).

Books, Proceedings

1. K. Kurokawa, M. Ishii, Y. Suda, A. Ichihara, and A. Nakano: “Live cell visualization of Golgi membrane dynamics by super-resolution confocal live imaging microscopy (SCLIM),” *Methods Cell Biol.* in press.
2. A. Nakano: “Super-resolution confocal live imaging microscopy – cutting-edge technology in cell biology,” *Proceedings of the IEEE Engineering in Medicine and Biology Society*, in press.
3. K. Hibino, M. Hiroshima, Y. Nakamura, and Y. Sako: “Single-molecule imaging measurements of protein-protein interactions in living cells,” An Integrated View of the Molecular Recognition and Toxinology – From Analytical Procedures to Biomedical Applications, Chapter 17, 432–453(2013).
4. M. Tominaga and A. Nakano: “Plant specific myosin XI, a molecular perspective,” *Frontiers Plant Cell Biol.* 3, 1-11 (2012).
5. C. Saito: “Ultrastructural analysis of single-membrane-bound organelles in plant cells: A comparison of high-pressure freezing versus chemical fixation,” *Plant Morphology*, 25, 11-14 (2013).
6. N. Imamoto, and T. Funakoshi: “Nuclear pore dynamics during the cell cycle.” *Curr. Opin. Cell Biol.* 24, 453-459 (2012). *Invited review* (査読有り)
7. N. Imamoto, and S. Kose: “Heat-shock stress activates a novel nuclear import pathway mediated by Hikeshi,” *Nucleus*, Sept 1; 3(5) [Epub ahead of print] (2012). *Invited review* (査読有り)
8. 廣島通夫, 佐甲靖志: “超解像顕微鏡 受容体イメージング” 現代化学 503, 56-57 (2013).
9. 佐甲靖志: “動的に生まれる膜 受容体の機能分担: 細胞内 1 分子計測技術で見たheregulin受容体の 2 量体形成反応” 日薬理誌 141, 240-244 (2013).
10. 山本明弘, 廣島通夫, 佐甲靖志: “Localization microscopyの原理と受容体イメージングへの応用” 光技術コンタクト 51, 33-40 (2013).
11. 小瀬真吾、今本尚子：“Hikeshi”は分子シャペロン Hsp70 を核に輸送し熱ストレスによる核の損傷から細胞を守る, 新着論文レビュー, 6月 12 日 (2012).
12. 小瀬真吾、今本尚子：新しい核-細胞質間運搬体分子“Hikeshi”的同定, 実験医学 30, 2621-2624 (2012).